

# Bracket AED 3

Wall Mount System for ZOLL AED 3™



- Low installation depth
- Lightweight
- Easy closing mechanism
- Two models Available
  - suitable for ZOLL AED 3™ Standard Carry Case
  - suitable for CRS medical Softbag AED 3
- 10G tested (dynamic)

### Technical Data for Standard Carry Case

- Dimensions (HxWxD) 351,5x258x183 mm
- Weight 3600 g

### Technical Data for CRS Softbag AED 3

- Dimensions (HxBxT) 334,5x251x153 mm
- Weight 3050 g



Bracket AED 3 with CRS Softbag AED 3



Bracket AED 3 for BLS Bag with ZOLL AED 3™ Standard Carry Case



CREATE SOLUTIONS - SAVE LIVES

**Bracket AED 3**  
for CRS medical Softbag AED 3

PN: 1042100003

**Bracket AED 3 for BLS Bag**  
for ZOLL AED 3™ Standard Carry Case

PN: 1042100004

**System AED 3**  
Bracket AED 3 and CRS Softbag AED 3

PN: 1041100003



\*Specifications subject to change without notice